



SANI-CLOTH® AF3

GERMICIDAL DISPOSABLE WIPE

Protecting your patients has never been easier. Fragrance Free—Ideal for use around patients and staff, especially those with respiratory sensitivities.

Sani-Cloth® AF3 Benefits:

- Perfect for equipment and other hard non-porous surfaces sensitive to alcohol
- Fragrance, alcohol, bleach, acid, phenol, acetone, and ammonia free
- Meets CDC, OSHA and CMS Tag F441 guidelines
- Bactericidal, Tuberculocidal, Virucidal



EPA Reg. No. 9480-9

Contact Time (Minutes)



of Microorganisms Effective Against

44

Disinfection Formulation

Alcohol Free Quat

Effective against 44 microorganisms in 3 minutes including the following MDROs, bloodborne pathogens and viruses:¹

- *Klebsiella pneumoniae* — NDM-1
- *Klebsiella pneumoniae* — Carbapenem Resistant
- ESBL producing *E. coli*
- CA-MRSA USA 300 and 400
- RSV
- HBV
- HCV
- HIV

SANI-CLOTH®
The #1 brand of disinfectant wipes in healthcare.²

¹See Technical Data Bulletin for a complete list of microorganisms. ²GHX Intelligence Data 2010.



ACCESSIBILITY: Deep well lid provides ample space to store and access next wipe.



COMPLIANCE: Prominent contact time for easy identification and improved compliance with accreditation standards.



IDENTIFICATION: Categorization of pathogens allows for quick identification of relevant efficacy claims.

ORDER INFORMATION

SANI-CLOTH® AF3	Reorder No.	Wipe Size	Packaging	Case Wgt. (lbs.)	Case Cube (ft.)	Tl/Hl
A. LARGE CANISTER	P13872	6" x 6.75"	160/can, 12 cans/cs	26.40	1.41	10/3
B. EXTRA-LARGE CANISTER	P63884	7.5" x 15"	65/can, 6 cans/cs	15.30	1.010	10/4
C. LARGE INDIVIDUAL PACKETS	H59200	5" x 8"	50/box, 10 boxes/cs	7.99	0.540	20/4
D. EXTRA-LARGE INDIVIDUAL PACKETS	U27500	11.5" x 11.75"	50/box, 3 boxes/cs	7.60	0.500	24/3
E. PORTABLE PACK	M928S80	8.2" x 9.8"	100/unit, 9 units/cs	22.63	0.948	10/4
F. PAIL	P1450P	7.5" x 15"	160/pail, 2 pails/cs	12.4	0.953	10/5
G. REFILL	P2450P	7.5" x 15"	160/refill, 2 refills/cs	11.1	0.529	18/4



SANI-CLOTH® EDUCATIONAL & COMPLIANCE ACCESSORIES

Developed to assist you in facility protocol compliance, as well as meet The Joint Commission and CMS requirements.^{3,4}



Infection Prevention Pak for mobile infection prevention – hold packets, gloves and other PPE



PDI 3-in-1 Universal Sani-Bracket® for wall/pole mounting of a Sani-Cloth® Canister



Sani-Canister Caddy® for placement on flat surfaces (bilingual)



Sani-TAG™ Equipment ID System for identification of clean/soiled equipment



Educational wall charts and signage (bilingual)

www.pdihc.com
1-800-999-6423

©2014 PDI, Inc. All rights reserved.

³Infection Prevention and Control (IC) Standard IC.02.02.01 and Elements of Performance (EP) 1 and 2, The Joint Commission, October 2009. ⁴Guideline for Disinfection and Sterilization in Healthcare Facilities, Centers for Disease Control and Prevention, 2008. Sani-Cloth® and Sani-Bracket® are registered trademarks of Professional Disposables International, Inc. Sani-Canister Caddy™ and Sani-TAG™ are trademarks of Professional Disposables International, Inc. PDI04125798

