


Victoria Bay

Quality products for everyday needs.

Med 1



DESCRIPTION

Med 1 is a blend of surfactants and solvents to lift stains like shoe polish, lipstick, and medicine. Med 1 is included in a four part kit designed to tackle all types of tough laundry stains. **NOT AVAILABLE IN ALL STATES.**

Please read entire container label for detailed instructions and safety information.

DIRECTIONS

1



Wear safety glasses & gloves
Usar gafas de seguridad y guantes

2



Pre-test fabrics in a hidden area
Probar el producto sobre un área de la tela que no esté a la vista

3



Spray liberally onto stain and agitate
Rocíe abundantemente sobre la mancha y agite

4



Wait 5 to 10 minutes
Esperar de 5 a 10 minutos

5



Laundry with normal procedures
Lave la ropa con los procedimientos normales

FEATURES

- Blend of surfactants

ADVANTAGES

- For use on: lipstick, shoe polish, and medical stains

BENEFITS

- Lifts stains

SPECIFICATIONS

Color Colorless
Fragrance Characteristic
pH Range 12.5 - 13.0
Dilution RTU
Density 8.67

REORDER

G00151 6/1 quarts

victoriabayproducts.com | 255 Route 1 & 9, Jersey City, NJ 07306 | (800) 535-5053

Guarantee: We guarantee your complete satisfaction with this product's performance. If dissatisfied for any reason, return any unused portion for credit within one year of the date of manufacture.